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
This document provides guidance on how to collect a fecal specimen and transfer to transport medium for diagnosis of acute bacterial diarrheal disease. Specimens that cannot be cultured within 2 hours of collection should be placed in Cary-Blair transport medium and refrigerated immediately.

Supplies needed:

- Specimen collection cup
- One tube of Cary Blair transport medium
- Sterile cotton or polyester-tipped applicators (swabs)
- Gloves
- Permanent marker
- Specimen labels or adhesive tape

1. Chill the tube of Cary Blair transport medium by placing it on ice packs or in the refrigerator 1 - 2 hours before collecting the specimen.

2. Gloves should be worn at all times when collecting and handling the specimen.

3. Collect stool from patients in a clean (no disinfectant or detergent residue) container with a tight-fitting, leak-proof lid. **Note:** Specimens should not be collected from bedpans, as they may contain residual disinfectant or other contaminants.

Alternatively if a stool specimen cannot be produced, rectal swabs may be collected as described in Step 5 below.

4. Remove the wrapper from the handle end of the sterile swab. Do not touch the tip of the swab.

5. **For Collected Stool Specimens:** Collect a small amount of stool by inserting a sterile cotton or polyester-tipped swab into the collected stool and rotating it. Mucus and shreds of intestinal epithelium, if present, should be sampled with the swab.

For Rectal Swabs: Moisten the swab in sterile transport medium, insert the swab through the rectal sphincter 2-3 cm (1-1.5 inches), rotate, and withdraw. Examine to ensure there is fecal material visible on the swab.

6. Immediately insert the swab into transport medium. The swab should be pushed completely to the bottom of the tube of transport medium.

7. Break off the top portion of the stick that was in contact with the gloved hand.

8. Repeat steps 4-7 for an additional sterile swab. Place the second swab in the same tube of Cary-Blair. Twist cap tightly on Cary-Blair tube to prevent leakage.

9. Record specimen number, and if possible the patient's name and date of collection, on the specimen label or adhesive tape and adhere to tube.

10. Safely dispose all contaminated materials. Do not reuse.

Doc. No. NERC/TEJ.002 Ver. No. 01 Effective Date: 05/01/2017

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

¹ Annex 2B- Case-based laboratory reprint from https://www.afro.who.int/sites/default/files/2017-Q6/IDSR-Technical-Guidelines_Final_2010_0.pdf

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
<https://www.cdc.gov/cholera/pdf/englishjobaid-stool-col-to-caryblair2.pdf>

Specimen Packaging and Domestic Transportation for Laboratory Confirmation of *Vibrio cholerae* 01/0139

FAECAL SPECIMENS CONDITIONING: 4 possible options			
Use gloves and lab coat when handling samples at all times. Specimen Label: Carefully identify specimens and indicate (using a permanent marker) patient name, date of collection, time, location of sampling and location of patient when likely infected. Lab Form: Use Annex 2B-105R case-based laboratory reporting form.*			
Faecal Specimen in Stool Container	APW (alkaline peptone water)	Wet and Dry Filter Paper (WFP/DFP)	Cary Blair medium, Faecal Sample or Rectal Swab
 Keep in initial stool container.	 Transfer faecal material from initial container into APW tube. NOTE: The faecal material should not exceed 10% of the volume of the APW enrichment.	 WET FILTER PAPER (WFP) Dip filter disk into wet fecal material with single-use device (forceps, needle), transfer into tube, add 2 to 3 drops of saline, close tube. DRY FILTER PAPER (DFP) Deposit one drop of wet fecal stool into filter paper. Air dry paper before placing into individual pouch with desiccant.	 For faecal samples: dip swab in liquid stool and transfer into Cary Blair medium. Rectal swab: Place swab directly into Cary Blair. No further manipulation is required.
Compatibility with testing methods (either directly from sample or after incubation steps in APW for those marked with *).			
RDT, culture, molecular analysis	RDT, culture, molecular analysis	WFP: culture, molecular analysis, RDT* DFP: molecular analysis	Culture, molecular analysis* and RDT*
MATERIAL REQUIRED			
Stool container (plastic, screw cap, 50ml, without disinfectant)	APW tubes with screw cap, transfer pipettes or swabs	WFP: filter paper discs (6mm Ø, non-sterile), saline solution, forceps or needle, 2ml tube (screw cap) DFP: Whatman cards (903 protein saver, FTA Elute Micro Cards), disposable transfer pipettes, individual pouches, desiccant	Cary Blair (semi-solid, bottles/tubes), swab (sterile, cotton/polyester)
Parafilm or sealing tape to seal packages and prevent leakage (not required for dry filter paper).			
CONSERVATION			
Ambient temperature (ideally 22-25°C). Do not refrigerate. Keep stool container out of direct sunlight.			
2 hours max. if delay > 2h, use Cary Blair.	Less than 24 hours	WFP: ideally less than 15 days DFP: no limitation	Follow manufacturer's instructions, on average 7 days
DOMESTIC TRANSPORTATION (national shipment, by road)			
Primary Containers  Primary containers are individually wrapped, and absorbent material placed between the primary containers and the secondary containers.	Absorbent Materials  Culture Media Cotton Balls Paper Towels	Samples are categorized "biological substances" category II. The use of triple packaging with UN3373 labels are required, alternatives are shown on the left. Samples must travel with corresponding documentation (lab request form and/or lab log). Include any results that may have already been performed, such as RDT results. Do not write the name of the organism on the outside of the package, only on the paperwork inside the box where appropriate. IMPORTANT: Indicate complete address and phone number for sender and recipient. Inform recipient laboratory about upcoming arrival of samples. TRANSPORT AT AMBIENT TEMPERATURE ICHC, April 2019 <small>View the faecal specimen packaging: https://www.gtfcc.org/sites/default/files/2019-04/04-2019-faecal-specimen-packaging_2019_04.pdf</small>	
Secondary Container  Sealed Polyethylene Container (1-inch-thick minimum)	Tertiary Containers  Plastic Bag Plastic Container Screw Cap Can Rigid container without top pack Corrugated Fibre Board UN 3373 Biological Substances Category II		

JOB AID: HOW TO COLLECT A FECAL SPECIMEN AND TRANSFER TO TRANSPORT MEDIUM

This provides guidance on how to collect a fecal specimen and transfer to transport medium for diagnosis of acute bacterial diarrheal disease. Specimens that cannot be cultured within 2 hours of collection should be placed in Cary-Blair transport medium and refrigerated immediately.

Supplies needed:

- One tube of Cary Blair transport medium
- Sterile cotton-tipped applicators (swabs)
- Gloves
- Permanent marker
- Specimen labels or adhesive tape

- Chill the tube of Cary Blair transport medium by placing it on ice packs or in the refrigerator 1 - 2 hours before collecting the specimen.
- Gloves should be worn at all times when collecting and handling the specimen.
- Collect stool from patients in a clean (no disinfectant or detergent residue) container with a tight-fitting, leak-proof lid. **Note:** Specimens should not be collected from bedpans, as they may contain residual disinfectant or other contaminants.
- Remove the wrapper from the handle end of the sterile swab. Do not touch the tip of the swab.
- Collect a small amount of stool by inserting a sterile cotton- or polyester-tipped swab into the stool and rotating it. Mucus and shreds of intestinal epithelium if present, should be sampled with the swab.
- Immediately insert the swab into transport medium. The swab should be pushed completely to the bottom of the tube of transport medium.
- Break off the top portion of the stick that was in contact with the gloved fingers.
- Repeat steps 5-7 for an additional sterile swab. Place the second swab in the SAME tube of Cary-Blair. Twist cap tight on Cary-Blair tube and specimen cup to prevent leakage.
- Adhere specimen label to the container or write on adhesive tape and secure to tube.
- Safely dispose all contaminated materials. Do not reuse.

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Antimicrobial Susceptibility Testing for Treatment and Control of Cholera

OBJECTIVE : To provide instruction for determining in vitro susceptibility of *Vibrio cholerae* 01/0139

METHODS

Combination of two methods:

- Agar disk diffusion method with antibiotic impregnated disks at predetermined concentrations
- Measurement of minimum inhibitory concentration (MIC) by using test-strips impregnated with a gradient of predefined concentrations of antibiotics

NOTE: Test strips are recommended for antibiotics for which no breakpoint is defined or when complementary tests are needed.
NOTE: Control strain(s) should always be set up in parallel with test strains.

MATERIALS REQUIRED

- Mueller Hinton Agar (MHA) plates (4 mm ± 0,5mm deep)
- Sterile saline solution (0,85% or 0,9%) ± test tubes of identical size for bacterial suspensions and the McFarland turbidity standard
- Sterile cotton tipped swabs
- Automatic disk dispenser or template with 5 or 6 disk spacing pattern and forceps
- Metric ruler (that can measure in mm)
- 0.5 McFarland turbidity standard
- Sheet of white paper with sharp black lines (can be prepared by hand or printed out)
- Control strain : *Escherichia coli* ATCC 25922

Antibiotic to be tested and recommended method:

- Azithromycin (AZ), MIC measurement
- Ciprofloxacin (CIP), MIC measurement (used to monitor the level of susceptibility to CIP for strains resistant to NA)
- Nalidixic acid (NA), disk diffusion (30 µg)* (used as an indicator of reduced fluoroquinolones susceptibility)
- Tetracycline (TE), disk diffusion (30 µg)


AZ, CIP and TE are the three selected antibiotics recommended for treatment of cholera according to GFCC: https://www.who.int/cholera/task_force/use-of-antibiotics-for-the-treatment-of-cholera.pdf?ua=1
• Store antibiotic disks and test-strips between -20°C and 8°C according to manufacturer's instructions.
• Check expiration date of antibiotic disks and test strips prior use.

PROCEDURE FOR DISK AND STRIP TESTING

1. Preparation of inoculum


Prepare a bacterial suspension with a few well-isolated colonies from an overnight (18-24 hours at 35 ± 2°C) agar culture in sterile saline solution adjusted to 0,5 McFarland by comparison to the standard.

NOTE: Ensure that the Standard is aliquoted into a tube that is the same size as the tube used to prepare the test suspension.

2. Inoculation of MHA

- Not more than 15 minutes after preparing the inoculum suspension.
- Dip cotton swab in bacterial suspension; remove excess liquid by pressing the swab against the inside wall of the tube.
- Streak the entire surface of the plate 3 times, rotating 60 degrees each time.
- Ensure the surface is completely dry before the next step.

3. Application of antibiotic disks

- Not more than 15 minutes after swabbing.
- Place the disks individually with an automatic disk dispenser or sterile forceps, gently pressing down onto the agar.
- Do not move disks once deposited.
- Replace lid, invert the plates and place in the incubator.


NOTE: Allow disks to reach ambient temperature before opening cartridge or container for storage.

4. Application of test-strips

- Not more than 15 minutes after swabbing.
- Place the strips on the agar according to the recommendations of the manufacturer.
- Do not move strips once deposited.
- Replace lid, invert the plates and place in the incubator.

NOTE: Test strips must be consistently stored in a freezer at -20°C. Allow strips that will be used to reach ambient temperature before placing on the agar.

5. Incubation: 18 hours at 35°C ± 2°C.

6. Reading: After 18 hours, observe the plate and measure the diameter (mm) of the inhibition ring. Read MIC value (in µg/ml) at the intersection of the lower part of the ellipse-shaped growth inhibition area with the test-strip. If a MIC value is between two fold dilutions, always round up to the highest value.

INTERPRETATION OF RESULTS

Results from the tetracycline disk should be used to predict susceptibility to doxycycline.
*NA susceptibility testing is an optional step but has proven to be a reliable and inexpensive first line screening tool for the detection of strains that should be individually tested for susceptibility to CIP. Strains categorized as sensitive to NA can be categorized as sensitive to CIP. If resistant to NA, the isolate should be tested for susceptibility to CIP by MIC measurement.

NOTE: Additional antibiotics can be tested for the epidemiological monitoring of strains (i.e., colistin, polymyxin B, ampicillin, chloramphenicol, trimethoprim/sulfamethoxazole, ceftazidime, first and second generations streptomycin).

Quality Control

If the control strain results are unexpected or out of range, any results on VC strains are invalid and the laboratory should investigate the sources of error.

INTERPRETATION: please refer to one of the following standards:

EUCAST: http://www.eucast.org/fileadmin/src/media/PDFs/EUCAST_files/Breakpoint_tables/v_9.0_Breakpoint_Tables.pdf (using interpretative criteria for the Enterobacteriaceae)

CLSI: https://clsi.org/media/1450/m45ed3_sample.pdf

NOTE: The standards are reviewed on a regular basis, please check you are using an up-dated version.

This document is intended for use by Reference Laboratories. Please keep update until testing is complete or in accordance with the laboratory sample retention policy.

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©GTFCC, V1.0, March 2021